

The Scalable Sleep System that Grows with Your Needs

- Scalable from Level IV to Full Level I PSG
- Multi-Night Recording for HST
- LCD Screen for Patient Set-Up, Signal Quality Check & Viewing Traces
- Bio-Sensor Auto Start
- RemLogic Software Included
- Hi-Quality Microphone (optional)
- TX Proxy Adds Wireless Interface for DC Devices

Embletta MPR

- Nasal Flow
- SpO₂
- Position
- Activity
- Microphone (optional)




Embletta MPR PG

- Abdominal effort
- Thoracic effort
- Nasal Flow
- Snore
- SpO₂
- Heart rate
- 1 x ExG
- Position
- Activity
- 1 x DC channel
- Microphone (optional)



ST Proxy

- 2 EEG
- 2 EOG
- Bipolar Chin EMG
- ECG
- 2 X Leg EMG



ST+ Proxy

- 6 EEG
- 2 EOG
- Tri-polar Chin EMG
- ECG
- 2 X Leg EMG

TX Proxy

- Differential Pressure
- 6 x DC (-5/+5V),
- Ambient Light
- Ability to record Flow/Pressure/Leak Data from various ResMed Therapy devices

