

SLEEP



Natus Sleep Supplies

Sensors • Cannulas • Ancillary Sleep Supplies

SENSORS

Respiratory Effort

Natus XactTrace® Respiratory Effort Systems

XactTrace is a comprehensive family of respiratory inductance plethysmography (RIP) effort belt systems, providing a range of sizes and solutions, in both Universal and Embla versions.

System options include:

- 1 Reusable (pp.52, 53)
- 2 Single Use Cut-to-Fit (pp.53, 54)
- 3 Single Use Pre-Sized (pp.54, 55)

XactTrace Universal Kit – Belts, Kits & Components

- For use with non-Embla PSG systems
- Reusable, one size fits all
- Integrated battery power; no cumbersome in-line modules

1 XactTrace Universal Belts & Belt Kits*					
Part Number	Type/Description	Length	Connector	Quantity	Compatibility
1421091	XactTrace Universal Kit – Complete (includes abdomen & thorax belts, abdomen & thorax snap sensors, two universal cables)	n/a	n/a	1 kit/pkg	Non-Embla systems that accept 1.5mm touchproof connectors
1421080	XactTrace Universal Abdomen Kit (belt, sensor, cable)	n/a	n/a	1 kit/pkg	
1421081	XactTrace Universal Thorax Kit (belt, sensor, cable)	n/a	n/a	1 kit/pkg	
1421040	XactTrace Adjustable Universal Belt	27-79" (65-200cm)	n/a	1/pkg	
XactTrace Universal Snap Sensors					
1421026	XactTrace Universal Snap Sensor - Abdomen	n/a	n/a	1/pkg	Non-Embla systems that accept 1.5mm touchproof connectors
1421027	XactTrace Universal Snap Sensor - Thorax	n/a	n/a	1/pkg	
XactTrace Universal Cable					
1451080	XactTrace Universal Cable	92.5" (235cm)	1.5mm touchproof	1/pkg	Non-Embla systems that accept 1.5mm touchproof connectors

*Also compatible with Embla SD32 systems

XactTrace Reusable Belts for Embla systems

- Plugs directly into Embla PSG systems; no additional modules required
- Low profile lock assures patient comfort
- Small, medium and large sizes available

1 XactTrace Reusable Belts				
Part Number	Belt Type/Description	Cable Length/Connector	Quantity	Compatibility
1421032	Small Abdomen, 22" (57cm) to 41" (104cm)	27.5" (70cm), yellow 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
1421033	Small Thorax, 22" (57cm) to 41" (104cm)	15.7" (40cm), blue 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
1421034	Medium Abdomen, 33" (83cm) to 56" (142cm)	27.5" (70cm), yellow 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
1421035	Medium Thorax, 33" (83cm) to 56" (142cm)	15.7" (40cm), blue 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium

Table continued on next page



SENSORS

Respiratory Effort

XactTrace Reusable Belts for Embla systems (continued from previous page)

1 XactTrace Reusable Belts				
Part Number	Belt Type/Description	Cable Length/Connector	Quantity	Compatibility
1421036	Large Abdomen, 47" (120cm) to 76" (194cm)	27.5" (70cm), yellow 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
1421037	Large Thorax, 47" (120cm) to 76" (194cm)	15.7" (40cm), blue 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
1421038	Medium Thorax, 33" (83cm) to 56" (142cm)	78" (200cm), blue 2-pin key	1/pkg	Embla S4000/S4500 systems
1421039	Medium Abdomen, 33" (83cm) to 56" (142cm)	78" (200cm), yellow 2-pin key	1/pkg	Embla S4000/S4500 systems
1421041	Small Thorax, 22" (57cm) to 41" (104cm)	78" (200cm), blue 2-pin key	1/pkg	Embla S4000/S4500 systems
1421042	Small Abdomen, 22" (57cm) to 41" (104cm)	78" (200cm), yellow 2-pin key	1/pkg	Embla S4000/S4500 systems
1421043	Large Thorax, 47" (120cm) to 76" (194cm)	78" (200cm), blue 2-pin key	1/pkg	Embla S4000/S4500 systems
1421044	Large Abdomen, 47" (120cm) to 76" (194cm)	78" (200cm), yellow 2-pin key	1/pkg	Embla S4000/S4500 systems



XactTrace Single Use Cut-to-Fit Universal Belts, Kits & Components

- Configurations for Non-Embla PSG systems (see cut-to-fit options for Embla PSG Systems on Page 54)
- Customized sizing for each patient; especially useful for infants and bariatric patients
- Complete kit includes: abdomen single use lock, thorax single use lock, abdomen reusable processor, thorax reusable processor, two rolls XactTrace single use cut-to-fit belt material
- Purchasable options: abdomen and thorax locks, cut-to-fit belt material, RIP processors and kits (lock, processor, cut-to-fit belt material)

2 XactTrace Single Use Cut-to-Fit Belt Rolls				
Part Number	Type	Roll Length	Quantity	Compatibility
1421000	XactTrace Single Use Cut-to-Fit Belt Material	65.5' (20m)	1 roll/pkg	Most PSG systems
XactTrace Single Use Cut-to-Fit Belt Kits *				
Part Number	Type/Description	Cable Length/Connector	Quantity	Compatibility
1421050	XactTrace Universal Single Use Cut-to-Fit Abdomen Kit (lock, processor & 1 roll of belt material)	n/a	1 kit/pkg	Non-Embla systems that accept 1.5mm touchproof connectors
1421051	XactTrace Universal Single Use Cut-to-Fit Thorax Kit (lock, processor & 1 roll of belt material)	n/a	1 kit/pkg	
1421052	XactTrace Universal Single Use Cut-to-Fit Belt Kit - Complete	n/a	1 kit/pkg	
Locks for XactTrace Single Use Cut-to-Fit Belts				
Part Number	Type/Description	Cable Length/Connector	Quantity	Compatibility
1421020	Abdomen	27.5" (70cm), yellow 2-pin key	1/pkg	Non-Embla systems that accept 1.5mm touchproof connectors. Requires 1421047 or 1421048 processor.
1421022	Thorax	15.7" (40cm), blue 2-pin key	1/pkg	
RIP Processors for XactTrace Single Use Cut-to-Fit Belts				
Part Number	Type/Description	Cable Length/Connector	Quantity	Compatibility
1421047	XactTrace RIP Processor Sensor, Thorax	78" (200cm) 1.5mm touchproof	1/pkg	Non-Embla systems that accept 1.5mm touchproof connectors
1421048	XactTrace RIP Processor Sensor, Abdomen	78" (200cm) 1.5mm touchproof	1/pkg	



*Also compatible with Embla SD32 system

SENSORS

Respiratory Effort

XactTrace Single Use Cut-to-Fit Belts & Components for Embla systems


- Configurations for Embla PSG systems only (see Universal cut-to-fit options on Page 53)
- Customized sizing for each patient; especially useful for infants and bariatric patients
- Purchasable options: abdomen and thorax locks, cut-to-fit belt material and cable extensions



1421020



1421022



2 XactTrace Single Use Cut-to-Fit Belt Rolls

Part Number	Type	Roll Length	Quantity	Compatibility
1421000	XactTrace Single Use Cut-to-Fit Belt Material	65.5' (20m)	1 roll/pkg	XactTrace Single Use Locks

Locks for XactTrace Single Use Cut-to-Fit Belts

Part Number	Type/Description	Cable Length/Connector	Quantity	Compatibility
1421020	Abdomen	27.5" (70cm), yellow 2-pin key	1/pkg	Embla N7000, Embla S4000/S4500, Embletta (PDS, 128, X100, MPR), Embla titanium, RIP processor
1421021	Abdomen	27.5" (70cm), yellow 2-pin key with built-in Embletta Gold Inline Inductor	1/pkg	Embletta Gold only
1421022	Thorax	15.7" (40cm), blue 2-pin key	1/pkg	Embla N7000, Embla S4000/S4500, Embletta (PDS, 128, X100, MPR), Embla titanium, RIP processor
1421023	Thorax	15.7" (40cm), blue 2-pin key with built-in Embletta Gold Inline Inductor	1/pkg	Embletta Gold only

Cable Extensions for XactTrace Single Use Cut-to-Fit Belts

Part Number	Type	Cable Length/Connector	Quantity	Compatibility
1451070	Abdomen	78" (200cm), yellow 2-pin key	1/pkg	Embla S4000/S4500 systems
1451072	Thorax	78" (200cm), blue 2-pin key	1/pkg	Embla S4000/S4500 systems



011521



011522



XactTrace Single Use Pre-Sized Universal Belts & Components

- Configurations for Non-Embla PSG Systems
- Same high quality RIP technology as XactTrace reusable belts
- No threading of material into the lock is required
- Available in four convenient standard sizes
- Pre-attached snap terminals for quick, easy set up (snap sensors required; sold separately)

3 XactTrace Single Use Pre-Sized Universal RIP Belts

Part Number	Type/Description	Quantity	Compatibility
011516	XactTrace Pre-Sized Single Use Belt, Pediatric	20/box (10 packs of 2)	XactTrace Pre-Sized Single Use Snap Sensors
011517	XactTrace Pre-Sized Single Use Belt, Small	20/box (10 packs of 2)	
011518	XactTrace Pre-Sized Single Use Belt, Medium	20/box (10 packs of 2)	
011520	XactTrace Pre-Sized Single Use Belt, Large	20/box (10 packs of 2)	

Snap Sensors and RIP Processors for XactTrace Single Use Pre-Sized Universal Belts

Part Number	Type/Description	Quantity	Compatibility
011521	XactTrace Snap Sensor, Abdomen, 26" (66cm) cable, yellow 2-pin key connector	1/pkg	XactTrace Pre-Sized Belts & RIP Processors
011522	XactTrace Snap Sensor, Thorax, 16" (41cm) cable, blue 2-pin key connector	1/pkg	
1421048	XactTrace RIP Processor Sensor, Abdomen, 78" (200cm) cable, 1.5mm female touchproof connectors	1/pkg	XactTrace RIP Processors & most PSG systems that accept 1.5mm female touchproof connectors
1421047	XactTrace RIP Processor Sensor, Thorax, 78" (200cm) cable, 1.5mm female touchproof connectors	1/pkg	

SENSORS

Respiratory Effort

XactTrace Single Use Pre-Sized Belts & Components for Embla systems

- For use with Embla systems only
- Same high quality RIP technology as XactTrace reusable belts
- No threading of material into the lock is required
- Available in four convenient standard sizes
- Pre-attached snap terminals for quick, easy set up (snap sensors required; sold separately)



3 XactTrace Single Use Pre-Sized Belts				
Part Number	Size	Connector	Quantity	Compatibility
011516	Pediatric	XactTrace Snap Sensors in this table	20/box (10 packs of 2)	Embla N7000, titanium, Embletta (all models), Embla S4000/S4500 systems
011517	Small	XactTrace Snap Sensors in this table	20/box (10 packs of 2)	
011518	Medium	XactTrace Snap Sensors in this table	20/box (10 packs of 2)	
011520	Large	XactTrace Snap Sensors in this table	20/box (10 packs of 2)	
Snap Sensors for XactTrace Single Use Pre-Sized Belts				
Part Number	Description	Cable Length/Connector	Quantity	Compatibility
011521	Abdomen	26" (66cm), yellow 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
011522	Thorax	16" (41cm), blue 2-pin key	1/pkg	Embla N7000, Embletta (all models), Embla titanium
013224	Abdomen	78" (200cm), yellow 2-pin key	1/pkg	Embla S4000/S4500 systems
013225	Thorax	78" (200cm), blue 2-pin key	1/pkg	Embla S4000/S4500 systems

*013382 Embletta Gold Inline Inductor required, one snap per sensor.

XactTrace Single Use Pre-Sized Belt Sizing

Choosing a Belt Size

- XactTrace Single Use Pre-Sized Belts come in four sizes: Pediatric, Small, Medium and Large
- Belt length can be adjusted to assure a good fit
- Belts should fit patients snugly without being uncomfortably tight
- Measure each patient's circumference and use the table below as a guideline

3 XactTrace Single Use Pre-Sized Belt Sizing	
Size	Fits Patient Circumference (Size Range)
Pediatric	11-22 in. (28-56 cm)
Small	16-35 in. (40-89 cm)
Medium	30-55 in. (76-140 cm)
Large	40-75 in. (102-190 cm)



DID YOU KNOW?

Humans are the only mammals that willingly delay sleep

SENSORS

Respiratory Effort

Respiratory Inductance Plethysmography Belts & Cables

- 5-color-coded sizes fitting a wide range of patients
- Soft, flexible, reusable belts are machine washable, sterilizable
- Simple snap connectors and interface cables
- Plug-and-play replacement RIP interface cables



Respiratory Inductance Plethysmography Belts

Part Number	Size	Colors	Quantity	Compatibility
4140090	Infant Small 15.7" (40cm)	White/Gray	2/pkg	Most systems that accept 1.5mm touchproof connectors & RIP Interface Cables listed below
4140089	Pediatric Medium 23.6" (60cm)	White/Red	2/pkg	
4140088	Adult Large 35.4" (90cm)	White/Blue	2/pkg	
4140087	Adult XL 47.2" (120cm)	White/Black	2/pkg	
4140114	Adult XXL 59" (150cm)	White/Yellow	2/pkg	

Respiratory Inductance Plethysmography Interface Cables

Part Number	Type	Lead Wire Length	Connector	Quantity	Compatibility
9101G*	Abdomen	78" (200cm)	1.5mm touchproof	1/pkg	Grass Comet-PLUS PSG systems
9102G*	Chest	78" (200cm)	1.5mm touchproof	1/pkg	Grass Comet-PLUS PSG systems
9131	Abdominal	8" (20cm)	Key	1/pkg	Grass AURA PSG & AURA PSG LITE
9132	Chest	8" (20cm)	Key	1/pkg	Grass AURA PSG & AURA PSG LITE
4140095	Abdomen	78" (200cm)	1.5mm touchproof	1/pkg	Most systems that accept 1.5mm touchproof connectors
4140096	Chest	78" (200cm)	1.5mm touchproof	1/pkg	Most systems that accept 1.5mm touchproof connectors

*Reduced output signal for Grass PSG systems



Adjustable Respiratory Inductance Plethysmography Belt Kit

- One size fits all (Adult)
- Reusable & easy to clean
- High quality, comfortable material
- Kit 019-481100 includes two 019-481200 belts, one 4140095 abdomen interface cable and one 4140096 chest interface cable

Part Number	Type	Belt Dimensions	Quantity	Compatibility
019-481100	Adjustable RIP Kit, Adult	23.6" (60cm) minimum to 55" (140cm) maximum	1 kit/pkg	Most systems that accept 1.5mm touchproof connectors
019-481200	Adjustable RIP Belt Only, Adult	23.6" (60cm) minimum to 55" (140cm) maximum	1/pkg	Most systems that accept 1.5mm touchproof connectors & RIP Interface Cables listed above

SENSORS

Piezo Respiratory Effort Belts

Piezo Crystal Effort Sensor – Double Buckle

- One Piezo crystal effort sensor – double buckle
- One large elastic double buckle belt

Part Number	Type	Belt Length	Lead Wire Length	Quantity	Compatibility
1421102	Chest or Abdomen	Large - 1" (2.5cm) x 67" (170.2cm)	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
1370G-KIT	Chest or Abdomen	Large - 1" (2.5cm) x 67" (170.2cm)	78" (200cm)	1/pkg	Reduced output signal for Grass PSG systems that accept 1.5mm touchproof connectors
1330A-KIT	Chest or Abdomen	Large - 1" (2.5cm) x 67" (170.2cm)	78" (200cm)	1/pkg	Key connector for Grass AURA PSG & AURA PSG LITE
1421101	Abdomen	Large - 1" (2.5cm) x 67" (170.2cm)	27.6" (70cm)	1/pkg	Yellow (abdomen) key for Embla
1421100	Chest	Large - 1" (2.5cm) x 67" (170.2cm)	15.7" (40cm)	1/pkg	Blue (chest) key for Embla



Large Double Buckle Replacement Belt

Part Number	Size	Compatibility
AIE071	Large - 1" (2.5cm) x 67" (170.2cm)	Replacement belt for these kits: 1421101, 1421102, 1421100, 1370G-KIT, 1330A-KIT



Piezo Crystal Effort Sensor Kit – Velcro® Tab

- One Piezo crystal effort sensor – double tab
- Two elastic tab belts; one small and one pediatric

Part Number	Type	Belt Lengths	Lead Wire Length	Quantity	Compatibility
ACE116	Chest or Abdomen	Small 1" (2.5cm) x 35" (88.2cm) Pediatric 1" (2.5cm) x 18" (45.7cm)	78" (200cm)	1 kit/pkg	Most systems that accept 1.5mm touchproof connectors.
1388G-KIT	Chest or Abdomen	Small 1" (2.5cm) x 35" (88.2cm) Pediatric 1" (2.5cm) x 18" (45.7cm)	78" (200cm)	1 kit/pkg	Reduced output signal for Grass PSG systems that accept 1.5mm touchproof connectors



Piezo Film/Crystal Effort Sensor Kit – Double Loop

- One Piezo crystal effort sensor – double loop
- Two elastic belts; one large and one extra large

Part Number	Type	Belt Lengths	Lead Wire Length	Quantity	Compatibility
019-479600	Chest or Abdomen, Piezo film	Large 1.5" (3.8cm) x 52" (132.1cm) XL 1.5" (3.8cm) x 75" (190.5cm)	78" (200cm)	1 kit/pkg	Most systems that accept 1.5mm touchproof connectors
019-479800	Chest or Abdomen, Piezo film	Large 1.5" (3.8cm) x 52" (132.1cm) XL 1.5" (3.8cm) x 75" (190.5cm)	78" (200cm)	1 kit/pkg	Reduced output signal for Grass PSG systems that accept 1.5mm touchproof connectors
ACE114	Chest or Abdomen, Piezo crystal	Large 1.5" (3.8cm) x 52" (132.1cm) XL 1.5" (3.8cm) x 75" (190.5cm)	78" (200cm)	1 kit/pkg	Most systems that accept 1.5mm touchproof connectors



019-479600 & 019-479800
ACE114
L and XL belts included with all three kits

Infant Piezo Beltless Effort Sensor

- Apply with adhesive tape

Part Number	Size	Lead Wire Length	Quantity	Compatibility
005732	Infant	7' (2.1m)	1/pkg	Most systems that accept 1.5mm touchproof connectors



SENSORS

Airflow

Thermistor Airflow Sensor

- Flexible, rubber coated probes for oral (one) and nasal (two)

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
1420431	Single Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Embla SD32
1420435	Single Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Embla systems with key connectors
1420430	Single Channel Oral/Nasal	Adult	23.5" (60cm)	1/pkg	Embla systems with key connectors
019-480400	Single Channel Oral/Nasal	Adult	35.5" (90cm)	1/pkg	Embla systems with key connectors
019-480600	Single Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors

Thermistor Airflow Sensor

- Comfortable design does not use nasal prongs

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
005890	Single Channel Oral/Nasal	Adult	7' (210cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors

Thermocouple Airflow Sensor

- Flexible, rubber coated probes for oral (one) and nasal (two)

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
1401G	Single Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
019-480200	Single Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
019-481000	Two Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
019-480800	Two Channel Oral & Nasal	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
019-480900	Two Channel Oral & Nasal	Pediatric	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors



SENSORS

Airflow

Thermistor Airflow Sensor for Cannula

- Works with any adult nasal cannulas
- Comfortable design does not use nasal prongs

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
005891	Single Channel Oral/Nasal	Adult	7' (210cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors



Thermocouple Airflow Sensor for Cannula

- Hooks to any adult nasal cannula
- Adjustable prongs for stable recordings and patient comfort

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
1401G-H	Single Channel Oral/Nasal	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
1421-H	Single Channel Oral/Nasal	Adult	48" (122cm)	1/pkg	Key - for Grass AURA PSG & AURA PSG LITE



Thermocouple Airflow Sensor for Cannula

- Attaches to any adult nasal cannula

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
1413G	Single Channel Oral	Adult	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors



Pediatric Thermocouple Airflow Sensor

- Stiff bendable prongs allow optimal patient comfort and signal quality

Part Number	Type	Size	Lead wire Length	Quantity	Compatibility
ACE131	Single Channel Oral/Nasal	Pediatric	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
1420406	Single Channel Oral/Nasal	Pediatric	78" (200cm)	1/pkg	Embla systems with key connectors



Infant Thermocouple Airflow Sensor

- Placed on the bridge of the infant's nose with the elements positioned under the nostrils

Part Number	Type	Size	Lead Wire Length	Quantity	Compatibility
ACE176	Single Channel Nasal	Infant-Preemies	71" (180cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors



SENSORS

Airflow



Disposable Embla® BreathSensor® Airflow Thermistor

- Saves time and money by eliminating the need for cleaning and sterilization

Part Number	Type	Size	Quantity	Compatibility
4140118	Single Channel Oral/Nasal	Adult	24/pkg	All BreathSensor Cables listed below
4140119	Single Channel Oral/Nasal	Children & Small Adults	24/pkg	
4140121	Single Channel Oral/Nasal	Infant	24/pkg	
4140120	Single Channel Oral/Nasal	Preemie	24/pkg	

Requires cable from table below



BreathSensor Cables

- Cable includes a battery pack

Part Number	Length	Connectors	Quantity	Compatibility
1420446	78" (200cm)	Key connector	1/pkg	Embla S4500, Embla N7000, Embletta Gold
4140069	78" (200cm)	1.5mm touchproof	1/pkg	Embla SD32; most systems that accept 1.5mm touchproof connectors



CannuTherm™ Integrated Thermal/Pressure Sensor

- Attaches to cannula; special port holds the thermal sensor
- Eliminates need to tape separate sensors together
- Complies with AASM recommendations for nasal pressure and nasal/oral temperature monitoring

Part Number	Type	Length	Quantity	Compatibility
105778	CannuTherm	78" (200cm) lead wire	1 sensor & 5 cannulas/pkg	Cannulas designed for CannuTherm, All systems that accept 1.5mm touchproof connectors
105750	Replacement adult nasal cannula	7' (213.4cm) cannula length	25/pkg	



ThermiSense® Sensor Kits

Kit includes:

- Oral/Nasal Thermistor
- Nasal cannula with holder
- Oral/Nasal cannula with holder
- Filter

Part Number	Type	Sensor Lead Wire Length	Cannula Length	Quantity	Compatibility
1420124	Adult	84" (213.4cm)	84" (213.4cm)	1/pkg	Embla systems with key connectors
1420129	Adult	30" (76.2cm)	84" (213.4cm)	1/pkg	
1420126	Pediatric	30" (76.2cm)	84" (213.4cm)	1/pkg	
1420127	Pediatric	84" (213.4cm)	84" (213.4cm)	1/pkg	

SENSORS

Snore Sensors

Piezo Snoring Sensor

- Records vibration in the neck using a round, flexible thermo-plastic disc

Part Number	Type	Lead Wire Length	Quantity	Compatibility
1420610	Snore, Piezoelectric	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
1420615	Snore, Piezoelectric	23.5" (60cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
1420605	Snore, Piezoelectric	78" (200cm)	1/pkg	Embla systems with key connectors
1420600	Snore, Piezoelectric	23.5" (60cm)	1/pkg	Embla systems with key connectors
019-479700	Snore, Piezoelectric	78" (200cm)	1/pkg	Most Xltek systems that accept 1.5mm touchproof connectors



Dynamic Microphone Snoring Sensor

- Converts snoring sounds into an analog signal
- Soft and lightweight; provides good seal against ambient noise

Part Number	Type	Diameter	Lead Wire Length	Quantity	Compatibility
1200G	Snore, Dynamic Microphone	0.9" (2.3cm)	76" (190cm)	1/pkg	Most systems the accept 1.5mm touchproof connectors
019-480300	Snore, Dynamic Microphone	0.6" (1.5cm)	76" (190cm)	1/pkg	



DID YOU KNOW?

The sound of snoring is the result of air turbulence vibrating the structures in the upper airway; these sounds can range from 50 – 100 decibels

SENSORS

Body Position Sensors



DC Body Position Sensor Kit

- Indicates the following positions: prone, supine, left, right or sitting/standing up
- Includes large, adjustable, Velcro buckle band, 1" (2.5cm) width


Part Number	Type	Lead Wire Length	Quantity	Compatibility
ACE119	DC Body Position, Voltage Type	78" (200cm)	1/pkg	Most systems that accept DC 1.5mm touchproof connectors
1440011	DC Body Position, Resistance Type	78" (200cm)	1/pkg	Most systems that accept a DC 1.5mm touchproof connectors
1574-KIT	DC Body Position, Voltage Type	12" (30cm)	1/pkg	Key connector - Grass AURA PSG & AURA PSG LITE
1440021	DC Body Position, Resistance Type	78" (200cm)	1/pkg	Key connector - Embla S4500
4140073	DC Body Position, Resistance Type	78" (200cm)	1/pkg	Key connector - Embla Sandman SD20
019-481300	DC Body Position, Voltage Type	78" (200cm)	1/pkg	1/8" Male Plug - Grass Comet, Comet-PLUS
019-479900	DC Body Position, Voltage Type	78" (200cm)	1/pkg	1/8" Female Plug - Embla Sandman SD32



AC Body Position Sensor Kit

- Indicates the following positions: prone, supine, left, right or sitting/standing up
- Includes large, adjustable, Velcro buckle band, 1" (2.5cm) width

Part Number	Type	Lead Wire Length	Quantity	Compatibility
1440010	AC Body Position	78" (200cm)	1/pkg	Most systems that accept 1.5mm touchproof connectors
019-406500	AC Body Position	78" (200cm)	1/pkg	



DID YOU KNOW?

Ducks at risk of attack by predators are able to balance the need for sleep and survival, keeping one half of the brain awake while the other slips into sleep mode


SENSORS

Limb Movement Sensor Kits

Limb Movement Sensor Kits

- Piezo film sensor technology
- Includes a pair of elastic ankle straps

Part Number	Type	Lead Wire Length	Quantity	Compatibility
1770G-KIT	Limb Movement	78" (200cm)	1 sensor & 2 bands/pkg	Most systems that accept 1.5mm touchproof connectors
1773-KIT	Limb Movement	67" (170cm)	1 sensor & 2 bands/pkg	Key - Grass AURA PSG & AURA PSG LITE
ACE146	Limb Movement	112" (285cm)	2 sensors & 2 bands/pkg	Most systems that accept 1.5mm touchproof connectors



Pressure Flow/Snore

Piezo Pressure Flow/Snore Sensors

- Isolated element provides immunity from environmental artifacts
- Strong, stable signal

Part Number	Sensor Type	Lead Wire Length	Cannulas	Quantity	Compatibility
4140094	AC Pressure / AC Snore	8" (20.3cm)	Sensor only, no cannulas included	1/pkg	Most systems that accept 1.5mm touchproof connectors
14838-KIT	AC Pressure / AC Snore	8" (20.3cm)	five 2' Adult Nasal	1 kit/pkg	Key - Grass AURA PSG & AURA PSG LITE
14832	DC Pressure / AC Snore	8" (20.3cm)	Sensor only, no cannulas included	1/pkg	Key - Grass AURA PSG & AURA PSG LITE



Multi-Press AC/DC Pressure Flow/Snore Sensor Kit

- AC flow/snore, low level output
- DC air pressure output provides true differential air pressure between the + and - pressure inputs
- Adjustable flow and snore sensitivities
- Kit includes Pressure Module, AC & DC wire-sets and AA batteries

Part Number	Sensor Type	Lead Wire Lengths	Connectors	Quantity	Compatibility
019-479500	AC/DC Pressure, AC Snore	AC Flow 9" (23cm), Snore 9" (23cm), DC Flow 78" (200cm)	AC Flow & Snore touchproof, DC Flow 1/8" Male Plug	1 kit/pkg	<ul style="list-style-type: none"> - Most systems that accept 1.5mm touchproof and 1/8" male plugs - Most Nicolet, Grass & Xitek systems are compatible - Embla & Schwarzer systems require adaptors



CPAP Pressure Tubing

- Tube with a Luer Lock connector to measure CPAP mask pressure

Part Number	Type	Length	Quantity	Compatibility
1420116	CPAP pressure tube	23.6" (60cm)	1/pkg	Embla N7000, Embla S4000/S4500, Embletta MPR, Embletta Gold. Most pressure flow/snore sensors & most CPAP masks
1420117	CPAP pressure tube	78.7" (200cm)	1/pkg	Embla N7000, Embla S4000/S4500, Embletta MPR, Embletta Gold. Most pressure flow/snore sensors & most CPAP masks



SENSORS

Oximetry



Nonin® Xpod®

- Clinically validated for challenging conditions such as motion and low perfusion
- Compatible with all Nonin PureLight® Sensors

Part Number	Type	Cable Length	Quantity	Compatibility
1430000	Nonin Xpod Oximeter	18" (45.7cm)	1/pkg	Embla N7000/S7000, S4000/S45000, titanium; Embletta Gold, PDS
011291	Nonin Xpod Oximeter	18" (45.7cm)	1/pkg	Embletta MPR
XPOD	Nonin Xpod Oximeter	80" (2m)	1/pkg	Grass Comet and AURA 32-channel EEG



Sensor Clip for Nonin Xpod

- Keeps Nonin Flex Sensor firmly connected to Xpod

Part Number	Type	Quantity	Compatibility
1431530	Xpod sensor clip	1/pkg	Nonin Xpod (rectangular)



Sensor Clip for Nonin Xpod

- Keeps Nonin Flex Sensor connected to Xpod

Part Number	Type	Quantity	Compatibility
013400	Xpod sensor clip	1/pkg	Nonin Xpod (tapered)



Nonin PureLight® Reusable FlexSensor & Disposable FlexiWrap® Adhesive Sensors

- Combines convenience of disposable sensor with economy of reusable sensor
- Available in multiple sizes
- Includes a starter pack of 25 disposable FlexiWraps

Part Number	Type	Cable Length	Quantity	Compatibility
1431000	FlexSensor, Adult	3' (1m)	1/pkg	Nonin oximeter
1431002	FlexSensor, Adult	9.8' (3m)	1/pkg	
1431011	FlexSensor, Infant	3' (1m)	1/pkg	
1431010	FlexSensor, Neonatal	3' (1m)	1/pkg	



Part Number	Type	Quantity	Compatibility
1431500	FlexiWrap, Adult	25/pkg	1431000 & 1431002 Flex Sensors
1431502	FlexiWrap, Infant	25/pkg	1431011 Flex Sensors
1431501	FlexiWrap, Neonatal	25/pkg	1431010 Flex Sensors

SENSORS

Oximetry

Nonin Reusable Pulse Oximetry Extension Cables

Part Number	Type	Cable Length	Quantity	Compatibility
1430500	Oximeter Sensor Extension Cable	3' (1m)	1/pkg	Nonin oximeter
222-490300	Oximeter Sensor Extension Cable	9.8' (3m)	1/pkg	
222-490400	Oximeter Sensor Extension Cable	19.9' (6m)	1/pkg	



Nonin PureLight® Reusable Soft Sensors

- Quick and easy application
- Comfortable fit for continuous monitoring
- Durable and easy to clean

Part Number	Type	Cable Length	Quantity	Compatibility
1431003	Soft Sensor, Small	3' (1m)	1/pkg	Nonin oximeter
105591	Soft Sensor, Small	9.8' (3m)	1/pkg	
1431004	Soft Sensor, Medium	3' (1m)	1/pkg	
105592	Soft Sensor, Medium	9.8' (3m)	1/pkg	
1431005	Soft Sensor, Large	3' (1m)	1/pkg	
105593	Soft Sensor, Large	9.8' (3m)	1/pkg	



Nonin Reusable Articulated Finger Clip Sensors

- For spot checking and short term monitoring
- Comfortable, self-aligning grip
- Durable design; easy to clean

Part Number	Type	Cable Length	Quantity	Compatibility
1431020	Finger Clip Sensor, Adult	3' (1m)	1/pkg	Nonin oximeter



Nonin PureLight® Disposable Sensors

- Soft, cushioning material provides comfort and secure placement for extended monitoring


Part Number	Type	Cable Length	Quantity	Compatibility
105720XL	Disposable Sensor, Adult	3' (1m)	24/pkg	Nonin oximeter



Masimo LNOP® Disposable Adhesive Sensors

- Disposable sensors limit cross-contamination risks and decrease site pressure
- Excellent adherence, increasing monitoring effectiveness in high-motion areas

Part Number	Type	Quantity	Compatibility
103848	LNOP Adhesive Sensor, Adult	20/pkg	Xitek PSG systems with 103847 cable
104038	LNOP Adhesive Sensor, Pediatric	20/pkg	
104557	LNOP Adhesive L-shaped Sensor, Infant	20/pkg	
104558	LNOP Adhesive L-shaped Sensor, Neonate	20/pkg	



SENSORS


Oximetry



Masimo LNCS® Disposable Adhesive Sensors with Cable Attached

- Disposable sensors limit cross-contamination risks and decrease site pressure
- Excellent adherence, increasing monitoring effectiveness in high-motion areas

Part Number	Type	Cable Length	Quantity	Compatibility
008781	LNCS adhesive sensor, Adult	18" (45.7cm)	20/pkg	Xltek PSG systems with 008776 cable
008782	LNCS adhesive sensor, Pediatric	18" (45.7cm)	20/pkg	
008783	LNCS adhesive sensor, Infant	18" (45.7cm)	20/pkg	



Masimo Y-I Multisite Reusable Sensors & Wraps

- For patients greater than 2.2lbs (1kg)

Part Number	Type	Cable Length	Quantity	Compatibility
104512	LNOP Y-I Reusable Sensor	3' (0.9m)	1/pkg	Xltek with 103847 cable
008786	LNCS Y-I Reusable Sensor	3' (0.9m)	1/pkg	Xltek with 008776 cable
104514	Disposable Standard Wrap	n/a	100/pkg	Y-I Reusable Sensor
104513	Disposable Petite Wrap	n/a	100/pkg	Y-I Reusable Sensor
104515	Disposable Foam Wrap	n/a	12/pkg	Y-I Reusable Sensor



Masimo Reusable Sensors

- LNCS DBI reusable adult soft sensor with 3ft cable

Part Number	Type	Cable Length	Quantity	Compatibility
008777	LNCS reusable sensor, Adult	3' (0.9m)	1/pkg	Xltek with 008776 cable



Masimo Reusable Finger Clip Sensors

- Adjustable hinge allows sensor to fit various finger sizes
- For short term spot checking

Part Number	Type	Cable Length	Quantity	Compatibility
104039	LNOP finger clip sensor, pediatric	3' (0.9m)	1/pkg	Xltek with 103847 cable
104040	LNOP finger clip sensor, adult	3' (0.9m)	1/pkg	Xltek with 103847 cable
008785	LNCS finger clip sensor, pediatric	3' (0.9m)	1/pkg	Xltek with 008776 cable
008784	LNCS finger clip sensor, adult	3' (0.9m)	1/pkg	Xltek with 008776 cable



Masimo Patient Cables

Part Number	Type	Cable Length	Quantity	Compatibility
103847	Masimo LNOP Patient Cable	12' (3.66m)	1/pkg	Xltek PSG systems
008776	Masimo LNCS Patient Cable	14' (4.27m)	1/pkg	Xltek PSG systems

SENSORS

Oximetry

Nellcor™ OxiMax™ Disposable Pulse Sensor

- For continuous oxygen saturation and pulse rate monitoring

Part Number	Type	Cable Length	Quantity	Compatibility
4140190	Disposable Sensor, Adult	36" (91.4cm)	24/pkg	Sandman SD32/SD20



Nellcor Adult Oxygen Saturation Sensor, Two Piece

- Combines reusable cable with sterile, disposable adhesive sensor
- Cost-effective alternative for patients with adequate peripheral pulses

Part Number	Type	Cable Length	Quantity	Compatibility
4140186	Disposable Sensor/reusable cable, Adult	36" (91.4cm)	1 cable & 24 sensors/pkg	Sandman SD32/SD20



Nellcor Adult-Neonatal Oxygen Saturation Sensor with Wraps

- Easy to use wrap-around design

Part Number	Type	Cable Length	Quantity	Compatibility
4140188	Reusable Sensor, Disposable Wraps, Adult/Neonatal	36" (91.4cm)	1 sensor & 50 wraps/pkg	Embla/Nellcor OxiMax



Nellcor Reusable Clip Sensor

- Convenient clip design assures quality and durability

Part Number	Type	Cable Length	Quantity	Compatibility
4140189	Reusable Clip Sensor, Adult, Finger	36" (91.4cm)	1/pkg	Embla/Nellcor OxiMax



Nellcor Sensor Extension Cable

Part Number	Cable Length	Quantity	Compatibility
4140184	10' (3m)	1/pkg	Embla/Nellcor OxiMax



CANNULAS



Adult Nasal/Oral Pressure Cannulas

- Semi-rigid wire in oral trunk allows precise placement
- Ideal length when pressure sensor is bedside or similar distance from patient
- Disposable

Part Number	Type	Supply Tube Length	Quantity	Compatibility
1420100	Nasal/Oral	1' (0.3m) Terminates in a male Luer Lock connector without filter	25/pkg	Most pressure flow/snore sensors with standard Luer Lock connectors
1420101	Nasal/Oral	7' (2.1m) Terminates in a male Luer Lock connector without filter	25/pkg	



Adult Nasal Pressure Cannulas

- Ideal length when pressure sensor is bedside or similar distance from patient
- Disposable

Part Number	Type	Size	Supply Tube Length	Quantity	Compatibility
1420002	Nasal	Adult	15" (38cm) Tubing terminates in male Luer Lock connector without filter	50/pkg	Most pressure flow/snore sensors with standard Luer Lock connectors
1420011	Nasal	Adult	47" (120cm) Tubing terminates in male Luer Lock connector without filter	50/pkg	
1420013	Nasal	Adult	71" (181cm) Tubing terminates in male Luer Lock connector without filter	50/pkg	



Infant Nasal Pressure Cannulas

- For newborns to three months old
- Ideal length when pressure sensor is located very close to the patient
- Disposable

Part Number	Type	Distance Between Nares	Supply Tube Length	Quantity	Compatibility
1420040	Nasal	0.3" (7.6mm)	3' (0.9m) Tubing terminates in male Luer Lock connector without filter	50/pkg	Most pressure flow/snore sensors with standard Luer Lock connectors



Intermediate Infant Nasal Pressure Cannulas

- For patients three months to one year of age
- Ideal length when pressure sensor is located very close to the patient
- Disposable

Part Number	Type	Distance Between Nares	Supply Tube Length	Quantity	Compatibility
1420030	Nasal	0.36" (9.1mm)	3' (0.9m) Tubing terminates in male Luer Lock connector without filter	50/pkg	Most pressure flow/snore sensors with standard Luer Lock connectors

CANNULAS

Pediatric Nasal Pressure Cannulas

- For patients one year of age and older
- Disposable

Part Number	Type	Distance Between Nares	Supply Tube Length	Quantity	Compatibility
1420020	Nasal	0.475" (1.21mm)	1' (0.3m) Tubing terminates in male Luer Lock connector without filter	50/pkg	Most pressure flow/snore sensors with standard Luer Lock connectors
1420050	Nasal	0.475" (1.21mm)	7' (2.1m) Tubing terminates in male Luer Lock connector without filter	50/pkg	



Nasal Pressure & Gas Sampling Cannulas

- Disposable

Part Number	Type	Size	Cannula Lengths	Quantity	Compatibility
005923	Nasal	Adult	One 7' (2.1m) Airflow Pressure tube & one 7' (2.1m) End Tidal CO ₂ tube. Tubing terminates in male Luer Lock connector without filter	25/pkg	Most pressure flow/snore sensors with standard Luer Lock connectors
005925	Nasal	Pediatric	One 7' (2.1m) Airflow Pressure tube & one 7' (2.1m) End Tidal CO ₂ tube. Tubing terminates in male Luer Lock connector without filter	25/pkg	



Adult ThermiSense® Nasal Pressure Cannulas with Filter

- Replacement cannula with holder for ThermiSense diagnostic thermistor temperature sensing system
- Disposable

Part Number	Type	Cannula Length	Quantity	Compatibility
1420110	Nasal	7' (2.1m) Tubing terminates in male Luer Lock connector with filter	25/pkg	Adult ThermiSense nasal/oral airflow sensors



Adult ThermiSense Nasal Pressure Cannulas

- Replacement cannula with holder for ThermiSense diagnostic thermistor temperature sensing system
- Disposable


Part Number	Type	Cannula Length	Quantity	Compatibility
1420111	Nasal	7' (2.1m) Tubing terminates in male Luer Lock connector without filter	25/pkg	Adult ThermiSense nasal/oral flow sensors
1420112	Nasal	2' (0.61m) Tubing terminates in male Luer Lock connector without filter	25/pkg	



Adult ThermiSense Nasal/Oral Pressure & Gas Sampling Cannulas


- Replacement cannula with holder for ThermiSense diagnostic thermistor temperature sensing system
- Disposable

Part Number	Type	Cannula Lengths	Quantity	Compatibility
1420114	Nasal/Oral	One 7' (2.1m) Airflow Pressure tube & one 7' (2.1m) End Tidal CO ₂ tube. Tubing terminates in male Luer Lock connector without filter	25/pkg	Adult ThermiSense nasal/oral airflow sensors



Enlarged view of sensor

CANNULAS



Pediatric ThermiSense Nasal Pressure & Gas Sampling Cannulas

- Replacement cannula with holder for ThermiSense diagnostic thermistor temperature sensing system
- Disposable

Part Number	Type	Cannula Lengths	Quantity	Compatibility
1420125	Nasal	One 7' (2.1m) Airflow Pressure tube & one 7' (2.1m) End Tidal CO ₂ tube. Tubing terminates in male Luer Lock connector without filter	25/pkg	Pediatric ThermiSense nasal/oral airflow sensors



Adult ThermiSense Nasal Pressure & Gas Sampling Cannulas

- Replacement cannula with holder for ThermiSense diagnostic thermistor temperature sensing system
- Disposable


Part Number	Type	Cannula Lengths	Quantity	Compatibility
1420113	Nasal	One 7' (2.1m) Airflow Pressure tube & one 7' (2.1m) End Tidal CO ₂ tube. Tubing terminates in male Luer Lock connector without filter	25/pkg	Adult ThermiSense nasal/oral airflow sensors



Disposable Cannula Filters

- Luer Lock connectors

Part Number	Type	Quantity	Compatibility
15809	Filter	1/pkg	Most cannulas & pressure flow/snore sensors



DID YOU KNOW?

Sea otters hold hands when they sleep so they don't drift away from each other

ANCILLARY SLEEP SUPPLIES

Electrodes

Disposable Small Electrodes with Cable Attached

- Pre-gelled Ag/AgCl discs, contact area 22mm x 28mm
- Comfortable foam backing

Part Number	Type	Lead Length	Quantity	Compatibility
1410514	Surface	39" (1m)	10/pkg	Most systems that accept 1.5mm touchproof connectors
1410515	Surface	59" (1.5m)	10/pkg	Most systems that accept 1.5mm touchproof connectors



Disposable Adhesive Center Snap Square Electrodes

- Solid gel Ag/AgCl electrodes, 1" x 1" cloth-backed contact area with robust adhesive
- Ideal for sleep studies

Part Number	Type	Quantity	Compatibility
019-427000	Square Snap	15/pkg	Most systems that accept 1.5mm touchproof connectors; requires reusable snap leads



Disposable Adhesive Offset Snap Round Electrodes

- Circular shape, 1.9" (48mm) diameter with offset snap
- Wet gel Ag/AgCl sensor material, foam backing

Part Number	Type	Quantity	Compatibility
1410505	Circular with offset snap	25/pkg	Most systems that accept 1.5mm touchproof connectors; requires reusable snap leads



Disposable Adhesive Offset Snap Small Electrodes

- 22mm x 28mm with offset snap
- Wet gel Ag/AgCl sensor material, foam backing

Part Number	Type	Quantity	Compatibility
1410510	Oval with offset snap	25/pkg	Most systems that accept 1.5mm touchproof connectors; requires reusable snap leads



Disposable Adhesive Center Snap Round Electrode

- Circular shape, 1-7/8" (48mm) diameter
- Solid gel Ag/AgCl sensor material with vinyl backing

Part Number	Type	Quantity	Compatibility
1410500	Circular with center snap	30/pkg	Most systems that accept 1.5mm touchproof connectors; requires reusable snap leads



See Page 16 for additional disposable snap electrodes

ANCILLARY SLEEP SUPPLIES

Snap Leads



Reusable Snap Lead Pairs

- Designed for use with disposable snap electrodes
- Orange lead wire, orange button snap, orange touchproof connector

Part Number	Type	Lead Length	Quantity	Compatibility
1411517	Snap Leads	39" (1m)	2/pkg	Most systems that accept 1.5mm touchproof connectors
1411524	Snap Leads	79" (2m)	2/pkg	



Reusable Snap Lead Pairs

- Designed for use with disposable snap electrodes
- Black lead wire, black button snap, black touchproof connector

Part Number	Type	Lead Length	Quantity	Compatibility
1411520	Snap Lead	60" (1.5m)	2/pkg	Most systems that accept 1.5mm touchproof connectors
1411535	Snap Leads	120" (3m)	2/pkg	



Reusable Snap Lead Pairs

- Designed for use with disposable snap electrodes
- Black lead wire, white button snap, white touchproof connector

Part Number	Type	Lead Length	Quantity	Compatibility
1411530	Snap Leads	120" (3m)	2/pkg	Most systems that accept 1.5mm touchproof connectors



Reusable Snap Lead, Ground

- Designed for use with disposable snap electrodes
- Green lead wire, green button snap, green touchproof connector

Part Number	Type	Lead Length	Quantity	Compatibility
1411516	Snap Leads	39" (1m)	1/pkg	Most systems that accept 1.5mm touchproof connectors



Short Reusable Snap Leads

- Designed for use with disposable snap electrodes

Part Number	Type	Lead Length / Color	Quantity	Compatibility
1411541	Snap Lead	24" (0.6m) / Red	1/pkg	Most systems that accept 1.5mm touchproof connectors
1411548	Snap Lead	24" (0.6m) / Green	1/pkg	
1411542	Snap Lead	29.5" (0.75m) / Black	1/pkg	

ANCILLARY SLEEP SUPPLIES

Snap Leads

Reusable Snap Leads

- Designed for use with disposable snap electrodes
- Black lead wire, black button snap, black touchproof connector

Part Number	Type	Lead Length	Quantity	Compatibility
F-SL48	Snap Lead	48" (1.2m)	3/pkg	Most systems that accept 1.5mm touchproof connectors
F-SL60	Snap Lead	60" (1.5m)	3/pkg	
F-SL72	Snap Lead	72" (1.8m)	3/pkg	
F-SL96	Snap Lead	96" (2.4m)	3/pkg	
F-SL120	Snap Lead	120" (3.0m)	3/pkg	



Reusable Snap Leads Breakout Cable

- Color-schemed for use with the Embla X100 Proxy
- Twisted red, yellow and black lead wires/button snaps breakout at 18" (0.46m); total lead length 28" (0.71m)
- Twisted green, orange and blue lead wires/button snaps breakout at 30" (0.76m); total lead length 40" (1m)
- Lead wires terminate in 1.5mm touchproof connectors

Part Number	Type	Lead Length	Quantity	Compatibility
1411546	Snap Lead Breakout Cable	28" (0.71m) / 40" (1m)	1/pkg	Embla X100 Proxy



Velcro Straps

Velcro Elastic Straps

- Provides patient comfort while holding components securely in place

Part Number	Description	Length	Quantity	Compatibility
1421234	Small Velcro Elastic Strap	30" (76cm)	1/pkg	Embla N7000, all Emblettas & most patient-worn PSG systems
1421236	Medium Velcro Elastic Strap	45" (114cm)	1/pkg	
1421238	Large Velcro Elastic Strap	72" (183cm)	1/pkg	
1421230	Velcro Elastic Strap Kit; one each Small, Medium, Large	One each of above	1 kit/pkg	



Cable Wraps

- Ideal for bundling cables and wires

Part Number	Type	Quantity	Compatibility
1450021	Velcro Cable Wraps	12/pkg	Most PSG systems

