

TSA II | NeuroSensory Analyzer



The Essential Tool in Sensory Nerve Evaluation

“QST of the thermal modalities is the only clinical test for quantitative assessment of small caliber sensory nerve fiber function, the primary transmitters of pain sensation”

Medoc
advanced medical systems


TSA II | NeuroSensory Analyzer

“Warm and cold perception threshold abnormalities may indicate the earliest evidence of neurological dysfunction”


- Demonstrate small-fiber abnormality, despite negative EMG (large-fiber) findings
- Identify early onset and progression of functional neurological deficit in patients with Low Back Pain
- Evaluate both nerve root compression and inflammation
- Document small-fiber allodynia in Complex Regional Pain Syndrome (CRPS)
- Evidence underlying Central Nervous System impairment in your difficult to diagnose patients
- Early detection of diabetic neuropathy of small-caliber nerve fibers
- Evaluate effectiveness of growth factors, aldose-reductase inhibitors and other new therapies
- Automatic comparison of test results to normative reference data
- Assess entrapment and toxic neuropathies
- Utilize test programs designed to assist in detection of malingers
- Portable, fast, easy-to-use clinical tool
- Versatile WINDOWS based software with selection of test paradigms for clinical, research and drug trial requirements
- Easy data export to Excel
- Outcome-based color graphic reports with automated narrative discussion of test results for payors and referral network
- FDA Cleared



Vibratory testing of the foot



Pre & Post analysis of thermal pain thresholds



Bilateral comparison report

“Sensory Nerve Conduction studies sample only the peripheral nerve system, whereas thermal and vibratory QST evaluates the entire sensory axis”

CE 0473



See instructions for use

Medoc Advanced Medical Systems, U.S. • Compass Medical Technologies, Inc
 1502 West Highway 54 – Suite 404, Durham, North Carolina 27707
 phone: 919.402.9600 • fax: 919.402.9607 • medoc@mindspring.com

Medoc Ltd. Advanced Medical Systems • World Headquarters
 1 Ha'Dekel St. PO.Box 423 Ramat Yishai, Israel 30095
 phone: +972-4-903.8800 • fax: +972-4-903.8808 • medoc@medoc-web.com

EU Authorized Representative CEpartner4U BV
 Esdoornlaan 13 3951 DB Maarn, The Netherlands, office@CEpartner4U.nl
 phone: +31-343-442.524 • mobile: +31-6-516-536.26 • fax: +31-343-442.162

www.medoc-web.com

