# Brite

# Completely wearable multi-channel fNIRS device



Truly portable and covering e.g. the prefrontal cortex, motor cortex, or visual cortex.



Can be combined with other techniques such as EEG and tDCS.



Measures oxy-, deoxy-, and total hemoglobin concentration changes.



Integrated 9-axis orientation sensor.



Easy analysis of your data with our superior analysis software; OxySoft.



Online connection via Bluetooth or offline recording with high sampling rate



The Brite 23 is optimized for frontal cortex and Brite 24 is optimized for hair-covered area.

# Main applications are found in:

- Brain oxygenation monitoring
- Sports science
- Functional studies
- Cerebral studies
- Hyperscanning

0

# Interested?

Contact us at askforinfo@artinis.com

www.artinis.com +31 481 350 980 Einsteinweg 17 6662 PW Elst The Netherlands

artinis

Page 1 of 2 Brite Leaflet - rev. 1902

# References to wireless fNIRS

Buchheit, M., Bishop, D., Haydar, B., Nakamura, F. Y. & Ahmaidi, S. Physiological Responses to Shuttle Repeated-Sprint Running. Int. J. Sports Med. 31, 402–409 (2010).

Maidan, I. et al. Changes in oxygenated hemoglobin link freezing of gait to frontal activation in patients with Parkinson disease: an fNIRS study of transient motor-cognitive failures. J. Neurol. 262, 899–908 (2015).

Maidan, I., Bernad-Elazari, H., Giladi, N., Hausdorff, J. M. & Mirelman, A. When is Higher Level Cognitive Control Needed for Locomotor Tasks Among Patients with Parkinson's Disease? Brain Topogr. 30, 531–538 (2017).

Maidan, I. et al. The Role of the Frontal Lobe in Complex Walking Among Patients With Parkinson's Disease and Healthy Older Adults: An fNIRS Study. Neurorehabil. Neural Repair 30, 963–971 (2016).

Mirelman, A. et al. Increased frontal brain activation during walking while dual tasking: an fNIRS study in healthy young adults. J. NeuroEngineering Rehabil. 11, 85 (2014).

Nieuwhof, F. et al. Measuring prefrontal cortical activity during dual task walking in patients with Parkinson's disease: feasibility of using a new portable fNIRS device. Pilot Feasibility Stud. 2, (2016).

Shadgan, B., Reid, W. D., Gharakhanlou, R., Stpublisher-ids, L. & Macnab, A. J. Wireless near-infrared spectroscopy of skeletal muscle oxygenation and hemodynamics during exercise and ischemia. Spectroscopy 23, 233–241 (2009).

# Technical specifications

#### TECHNOLOGY

MEASURES DATA ANALYSIS SOFTWARE OPERATING SYSTEM LIGHT SOURCE CHANNELS

WAVELENGTHS DETECTORS OPTODE DISTANCE POWER

TOTAL WEIGHT SIZE

ENVIRONMENT ELECTROMAGNETIC COMPATIBILITY OPTIONAL OFFLINE DATA STORAGE RESOLUTION SAMPLE RATE ORIENTATION SENSOR NIRS + OTHER MODALITIES Continuous wave near infrared spectroscopy using the modified Lambert-Beer law Oxy-, deoxy-, and total hemoglobin OxySoft Windows 10 10 or 11 light emitting diodes, each with 2 wavelengths 23, 24, or up to 54 channels with dual Brite. Different setup's and short channels possible Standard nominal 760 and 850 nm, others possible 7x or 8x photodiodes with ambient light protection 30 mm Fast charging battery that lasts up to 3 hours on a single charge. Charging with powerbank possible. Approximately 300 grams including battery Battery housing: 75x75x30 mm. Headbands and head caps available in multiple sizes: kids version (3months - 6 years) and adults (XS - XL) Operating temperature: 10 - 35 °C No interference with EEG, ECG or EMG Can be combined with the PortaSync 100+ hours 16 bit 50 or 100 Hz 9-axis accelerometer We deliver the following packages: Brite + TMSi EMG package (2 or more channels) Brite + TMSi EEG package (16 channels or more) Brite + tDCS (STARSTIM)







## Get a quote:

## askforinfo@artinis.com

# What's in the box?

Brite research package

Brite 23 or 24 Strong & sturdy Pelicase

Battery charger Laptop User guide Bluetooth dongle OxySoft, data analysis software Neoprene headband/head cap

©2019 Artinis Medical Systems all rights reserved. All specifications and appearances can change without further notice. All other trademarks or registered trademarks are the property of their respective owners. Note: the Brite is used as a research device only.

www.artinis.com

Page 2 of 2 Brite Leaflet - rev. 1902