



## **ENT Dual Needle Electrode**

**Friendship Disposable ENT Dual Needle Electrodes** are a specially designed version of our Monopolar Needle Electrodes. They are primarily used during IOM for monitoring EMG activity during ENT surgery. The needles are color coded and come in boxes of 24 with six of each of the following colors: red, blue, orange or purple.

- · Individually packaged
- · Available in boxes of 24
- Sterile
- Single patient use only

## **Disposable ENT Dual Needle Electrodes**

C	Order No	Part No	Needle Diameter	Gauge	Needle Length	Hub Color	Recording Area	Packing pr Box
1	202	NE-D1-1500/12/0.4	0.4 mm	12 mm	150 cm	R,B,O,P	1.5 mm TP	24 pcs

<sup>\*</sup> Each Box contain 24 Individual packed sets in 4 colors;

6 pcs Red/Red twisted, 6 pcs Blue/Blue twisted, 6 pcs Orange/Orange twisted and 6 pcs Purble/Purble twisted.

