


**PATENT
PENDING**
EP 20214073.7
US 63/125,594



Introducing:
The Flow Arm
Quick coil placement and fixation

- Optimized for clinical practice
- Smooth, effortless movement
- Balanced coil weight
- No strain on the operator
- Suitable for all treatment coils



4 key features of the Flow Arm

3

Stepper

Allows a wide range of coil positions.

1

Angulator

Easy and accurate coil placement.

2

Balancer

Ensures contact between the coil and the patient's head and reduces the strain on the operator.

4

Slider

Allows smooth coil movement and improves patient flow.



The Flow Arm is to be mounted on the back of the treatment chair.
Older chairs may need an upgrade.

*Indications for use

MagVenture TMS Therapy® is FDA cleared for:

“Treatment of Major Depressive Disorder in adult patients who have failed to receive satisfactory improvement from prior antidepressant medication in the current episode and as an adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD).”

MagVenture TMS Therapy® is CE approved for:

“Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.”

www.magventure.com

MagVenture TMS Therapy
Stimulating mental health.