


The basic solution for repetitive Transcranial Magnetic Stimulation

MagPro R20



- The optimal choice when running a limited number of TMS sessions per day
- CE approved and FDA cleared for depression treatment*
- Quick and simple protocol set-up
- Up to 5,500 pulses in one session
- Mobile



Need more options?

The MagPro R20 Wi-Fi package contains a secure Wi-Fi connection and a configurator, providing you with more options:

- Protocol set up on an external computer.
- Wireless upload of the protocol and settings on the MagPro R20.
- Print documentation of performed sessions with your logo, header and patient details.



Technical data

Output power roll-off		Environmental Data		Power Supply	
Max. output at 5 pps	100 %	Operating temperature	10-30°C (50-86°F)	Mains voltage	230V~, 50/60 Hz. according to IEC 60601-1
Max. output at 10 pps	76 %	Storage temperature	5-50°C (41-122°F)	Max. power consumption	800 VA, operation from 100V~ or 120V~ through transformer.
Max. output at 15 pps	61 %	Operating humidity	30-60 % RH		
Max. output at 20 pps	54 %	Storage humidity	20-80 % RH		
Mechanical Data		Technical specifications		Ordering Number	
HxWxD	150 x 390 x 440mm 5,9 x 15,3 x 17,3 in	Maximum Repetition Rate	20 pps	9016E0861	
Weight	20kg (44 lbs)	Pulse mode	Standard		
		Waveform	Biphasic - 280 ms		

*Indication for use

MagVenture TMS Therapy® is FDA cleared for:

“Treatment of Major Depressive Disorder in adult patients who have failed to receive satisfactory improvement from prior antidepressant medication in the current episode.”

MagVita TMS Therapy® is CE approved for:

“Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.”

More info at

www.magventure.com

MagVenture
Research. Treatment. Results.