



# THE ADVANCED SOLUTION FOR TRANSCRANIAL MAGNETIC STIMULATION


MagPro X100



- ✓ For demanding repetitive protocols and high number of pulses
- ✓ Monophasic, biphasic, and half-sine waveforms
- ✓ Twin, dual & theta burst pulses
- ✓ Enhanced stimulation capabilities with MagOption
- ✓ CE approved for depression treatment\*

# NEED MORE OPTIONS?

- Data storage and reporting capability
- Real-time information on important performance parameters such as remaining stimuli and actual power delivered
- Complete packages available for interleaved TMS/fMRI
- Real double blinded research studies option



## Technical data



	MagPro X100	MagOption
<b>Mechanical Data</b>		
HxWxD	210 x 530 x 400mm 8.3 x 20.9 x 15.7 in	130 x 530 x 400mm 5.1 x 20.9 x 15.7 in
Weight	36 kg (79 lbs)	28 kg (62 lbs)
<b>Technical specifications</b>		
Maximum Repetition Rate	100 pps, 250 pps	100 pps, 250 pps
Pulse mode	Standard	Standard, Dual, Twin and Power
Waveform	Monophasic, Biphasic and Biphasic Burst (Theta Burst)	Monophasic, Biphasic, Biphasic Burst and Halfsine
Current direction	Normal or Reverse	Normal or Reverse
<b>Ordering Number</b>		
	9016E0711	9016E0731



**\*MagVenture TMS Therapy® is CE approved for:**  
 Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.

For subjects responding to the treatment of MDD, who have symptoms of comorbid anxiety, the clinical benefit can additionally include an improvement of anxiety symptoms.

Treatment of psychoactive substance use disorder (PSUD) in adult patients.

An adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD).

[www.magventure.com](http://www.magventure.com)

