



CEFALY

Migraines?

Cefaly® can help you!



**ACUTE AND PREVENTIVE
MIGRAINE TREATMENTS**

www.cefaly.com



What is Cefaly®?

Cefaly® is an External Trigeminal Nerve Stimulation device (e-TNS) for migraine treatment. A self-adhesive electrode is placed on the forehead and the Cefaly® device is magnetically connected to this electrode. Precise micro impulses are then sent through the electrode to the upper branch of the trigeminal nerve either to relieve the headache pain during a migraine attack (Acute setting) or to prevent future migraine attacks (Prevent setting).

Who is Cefaly® for?

Cefaly® is intended for patients who suffer from migraine with or without aura particularly when drug consumption needs to be reduced.

The **Acute setting** is indicated for the acute treatment of migraine attacks. The **Prevent setting** is indicated for migraineurs who have frequent attacks (3 or more per month) and need a migraine prophylaxis to reduce the frequency of attacks.

How is Cefaly® used?

To use Cefaly®, place a self-adhesive electrode on the forehead, then position the magnetic device on the electrode. The device is powered by a rechargeable battery and is controlled by means of a button on the front of the Cefaly®.

Acute setting should be used for a single 1-hour session during the migraine to stop the attack or relieve the headache pain.

Prevent setting should be used once daily for 20 minutes (preferably in the evening) to reduce the frequency of migraines attacks.




Cefaly® works!

The efficacy and safety of Cefaly® have been demonstrated by several clinical trials in the USA and Europe.

Acute treatment of migraine :

After a 1-hour session with the Acute setting, 32% of patients were pain free and 85% had pain relief.


% OF PATIENTS WHO
HAVE PAIN RELIEF



Prophylactic treatment of migraine :

After 3 months of daily 20-minute sessions with the Prevent setting, 3 out of 4 patients experienced a 50% or more reduction in migraine days and took on average 1/4 the amount of acute medication.

% OF RESPONDERS



Cefaly® is safe!

No serious adverse events were reported in clinical studies. Also, there have been no incidents or significant side effects in any clinical studies or in the tens of thousands of patients who have used Cefaly® for several years.

Test Cefaly®: Ask your doctor

For **CHF 75.–** (incl. VAT and shipping) you can test Cefaly for **40 days**.

You buy the Cefaly at the normal retail price of **CHF 421.–** (incl. VAT and shipping). After the trial period you have the following options:

1. You keep the Cefaly.
2. You return the Cefaly to us and we will refund you **CHF 346.–** (sales price CHF 421.– minus CHF 75.– rental fees).
3. You extend the test phase under

Cefaly® kit for a 40-day trial rental:

Brand new Cefaly® device

Storage case

Detailed instruction manual

Power adapter + cable

Cefaly® electrodes



Do you have any question? We answer you with pleasure. **Exclusive partner in Switzerland.**

NEUROLITE

Advanced Medical Solutions

Cefaly Switzerland c/o Neurolite AG

Hühnerhubelstr. 79
3123 Belp
Switzerland

+41 31 529 28 32
info@cefaly.ch
www.cefaly.ch

Facebook:
Cefaly Schweiz