



Corkscrew Subdermal Needle Electrode

Friendship Corkscrew Subdermal Needle Electrodes are a special version of the normal Subdermal Needle Electrode used to collect bioelectric signals (EEG, ECG and EMG) from the body during clinical examinations, treatments and surgical operations

The corkscrew design is particularly useful during anesthesia as it attaches like a small screw needle to the top layer of the scalp so it stays securely in place. Two varieties of the corkscrew needle electrode are available.



Corkscrew Subdermal Needle Electrode – Single

Friendship Disposable Corkscrew Subdermal Needle Electrodes – Single are a specially designed version of our Subdermal Needle Electrode intended for measuring EEG and EP signals, typically during interoperative monitoring. The corkscrew form is designed so it stays securely in place during surgery and examinations for clear and reliable results.

- Sterile
- Single use only
- Individually packaged

Corkscrew Subdermal Needle – Single

Item No.	Part No	Needle Diameter	Needle Length	Lead Wire Length	Lead Wire color	Connector	Вох
1001	NES-D1-1200/0.6	0.60 mm	N/A	120 cm	6 color	1.5 mm TP	24 pcs
1001-1	NES-D1-1200/0.6	0.60 mm	N/A	120 cm	12 color	1.5 mm TP	24 pcs
1003	NES-D1-1500/0.6	0.60 mm	N/A	150 cm	6 color	1.5 mm TP	24 pcs
1004	NES-D1-2000/0.6	0.60 mm	N/A	200 cm	6 color	1.5 mm TP	24 pcs
1005	NES-D1-2500/0.6	0.60 mm	N/A	250 cm	6 color	1.5 mm TP	24 pcs



Corkscrew Subdermal Needle Electrode

Friendship Europe ApS // Industrivej 3 // 4000 Roskilde, DENMARK // info@friendshipeurope.com // www.friendshipeurope.com





Corkscrew Subdermal Needle Electrode – Twisted Pair

Friendship Disposable Corkscrew Subdermal Needle Electrodes Twisted Pair is a specially designed version of our Subdermal Needle Electrode intended for measuring EEG and EP signals typically during interoperative monitoring. It has been developed to reduce possible noise issues during monitoring. The corkscrew form is designed to stay securely in place during surgery and examinations for clear and reliable results.

- Sterile
- Single use only
- Individually packaged

Corkscrew Subdermal Needles - Twisted

Item No.	Part No	Needle Diameter	Needle Length	Lead Wire Length	Lead Wire color	Connector	Вох
1006	NES-D2-1500/0.6	0.60 mm	N/A	150 cm	12 color	1.5 mm TP	12 Pair
1007	NES-D2-2000/0.6	0.60 mm	N/A	200 cm	12 color	1.5 mm TP	12 Pair

Extension Cables for 1,5 mm TP Lead Wire

Order No	Part No	Cable Length	Cable Color	EQ side	Electrode side	Packing pr Box
1060	EW1000	100 cm	P1 or P2	1.5 mm TP	1.5 mm femal	6 pcs
1061	EW1500	150 cm	P1 or P2	1.5 mm TP	1.5 mm femal	6 pcs
1062	EW2000	200 cm	P1 or P2	1.5 mm TP	1.5 mm femal	6 pcs

* Color set P1: (Red, Orange, Yellow, Green, Blue, Purple) Color set P2: (Black, White, Light Blue, Brown, Pink, Grey)



Corkscrew Subdermal Needle Electrode